HUMITAL SEQUIRON BIO

Mixture of trace elements boron (B), copper (Cu), iron (Fe), Manganese (Mn), molybdenum (Mo) and zinc (Zn)

IN THE FIELDS CHOOSE THE LIFE

FERTILIZERS Organics, Liquids, Correctives, Idrosolubles, Organ-Min

F.ne Rocchetta, 63 62027 San Severino Marche (MC) tel s10 0731355125 email: info@bhumeseril

CE FERTILIZER

Suitable for use in horticulture

DESCRIPTION

It is a formulated to based on chelated microelements, including the magnesium that is in a quantity high, made to afford all the microelements in the form readily assimilable.

This formulation is not a mixture of individual compounds with salt of magnesium, but it is complex reaction of chelation in the liquid phase and subsequently concentrated, dried and spraizzata to make the powder into a form that is highly water soluble.

Then inside of each microgranules are in perfect balance all the elements declared in the composition.

Thanks to a high presence of EDTA chelated microelements that ensure maximum absorption and do not create antagonism between them, it improves the equilibrium state of the plant and the organoleptic characteristics of the fruit (flavor, size).

METHOD 'AND DOSAGE

Typically is apply for fertigation at doses of 2-5 kg / ha.

Suitable for all crops VEGETABLES, TREES AND INDUSTRIAL FLORAL and cultures, and ORNAMENTAL NURSERIES,

repeating more times, in function on the cultivation needs.

The doses given are the result of applied and calculated for distributions to normal volumes of water. For the correct application in the specific climate, soil and crop, it is advisable to consult your service technician. Respect the time of shortage.

COMPOSITION

Boron (B) soluble in water	0,6%
Copper (Cu.) soluble in water	0.5%
Copper (Cu.) chelated with EDTA	0.5%
Iron (Fe.) soluble in water	4.0%
Iron (Fe.) chelated with EDDHA	2.7%
Iron (Fe.) chelated with EDTA	1.3%
Manganese (Mn.) soluble in water	1.0%
Manganese (Mn.) chelated with EDTA	1.0%
Molybdenum (Mo.) soluble in water	0.03%
Zinc (Zn.) soluble in water	0.8%
Zinc (Zn.) chelated with EDTA	0.8%

To be used only in the case of a recognized need. Do not exceed the appropriate dose

Copper-Manganese-Zinc: chelating agent EDTAIRON: EDDHA chelating agent EDTA -

Copper-Manganese-Zinc: chelating agent EDTA PH range that ensures a good stability of the chelated fraction: 3-9IRON: EDDHA chelating agent EDTA -

The methods of analysis not reported are internal methods of the manufacturer.



Choncimer srl - F.ne Rocchetta, 63 - San Severino Marche (MC) tel. 0733636325 : e-mail:info@choncimer.it : P.IVA 01265640431